


APS Adhesion Microscope Slide

Matsunami APS hydrophobic adhesion slide offers superior tissue adhesion, staining performance. Proprietary slide coating ensures consistent performance even after exposure to heat stress and humid environmental conditions during IHC staining.

Recommended for use on automated IHC platforms.

Matsunami microscope slides are manufactured in Japan using the high transparency glass, which is optically clear, exceptionally flat, low fluorescence.

Slides are packaged 100 per box, 10 boxes per case (1000 slides), 4 cases per carton.



Product No.	Product Name	Color	Size (mm)	Thickness (mm)	Frosted area width (mm)	Package	Remarks
SUAPS19	APS White	White	75×25	0.9 to 1.2	20×25	100 pcs	90 deg. corners

Result of comparative performance evaluation of APS costed slide.

IHC Platform: Benchmark Ultra

Antibody: Vimentin V9

Tissue section: Human mammary gland

Conventional hydrophobic slide	APS hydrophobic adhesion slide
